

Concentrated Liquid Chromogenic Culture Medium



Concentrated liquid chromogenic culture medium ST-SSP

Keep 5ml of the sample solution in ST-25, and then inject 1.65 ml of the concentrated liquid chromogenic medium (ST-SSP) completely into ST-25 for culture.

Performance Characteristics

ST-SSP can simultaneously detect coliforms and Escherichia coli for qualitative identification.

Name	Description	Packing	Validity period
ST-SSP	Coliforms and Escherichia coli for qualitative identification	10pcs/bag, 300pcs/box	1year

Operation and Coloring



Add SSP concentrated medium



Primary color – negative



Blue – coliforms positive



Fluorescence -- Escherichia coli positive