



Sterile Sampling Bag




Sterile Sampling Bag

Model	Dimension (mm)	Capacity (ml)	Material	Operating temperature	Packing/carton (pcs)	Sterilization Method	Validity Period
AL1276	127 X 76	60	Polyethylene	-50~85℃	500	Ethylene oxide sterilization	3 years
AL1776	178 X 76	120			500		
AL2211	229 X 114	530			500		
AL2214	229 X 140	700			500		
AL3018	300 X 180	1620			250		
AL2211B	229 X 114	530			500		
AL2214B	229 X 140	700			500		
AL3018B	300 X 180	1620			250		
AL2525B	254 X 254	1620			250		
AL3525B	356 X 254	2700			250		
AL5038B	508 X 380	5441			200		

Precautions:

1. Do not use for other application than sample sampling, collection, and storage.
2. Keep the sample temperature within the operating temperature range.
3. Used to freeze liquids after which handle and contact carefully.
4. When placing the sample directly in homogenization, be sure to release the sealing strip before homogenizing.
5. After use, remove the contents and dispose of them as specified.
6. In order to improve product quality, specifications and price changes will be notified separately.

A type and B type (with B letter) difference



Operation Instruction (B type)





Step 1



Step 2



Step 3



Step 4



Step 5



Step 6 Finished



B type sampling-finished effect



A type sampling-finished effect