


Swab Sampler



Performance characteristics:

- Easy to peel and pack;
- Quick leak-proof cover;
- High quality clean polyester pointed cotton swab;
- The cotton swab is firmly fixed to the cover;
- Validity period: 18 months.



Swab Sampler

Model	Specification	Application	Sterilization method	Packing	Validity period
ST-25PBS	Phosphate buffered water 10ml	Conventional strain sampling	X-ray irradiation sterilization	400pcs /carton	2years
ST-25BPW	1% buffered peptone water 10ml	Salmonella enrichment, the repair of oligotrophic bacteria and injured bacteria			
ST-26PBS	Phosphate buffered saline 10ml, extension type	Conventional strain sampling			
ST-26BPW	1% buffered peptone water 10ml, extension type	Salmonella enrichment, the repair of oligotrophic bacteria and injured bacteria			
ST-25MRD	0.1% phosphate buffered peptone saline	The repair of oligotrophic bacteria and injured bacteria			
AL-SWAB(P)	1% buffered peptone water 10ml	Salmonella enrichment, the repair of oligotrophic bacteria and injured bacteria			



ST-25LB	Lethen broth buffer 4ML/10ML	Enrichment culture of bacterial used in cosmetics or other substances (containing quaternary ammonium compounds or cationic surfactants)			
ST-26LB	Lethen broth buffer 4ML/10ML	Enrichment culture of bacterial used in cosmetics or other substances (containing quaternary ammonium compounds or cationic surfactants)			18 个月 18months

Note: Various buffers can be customized according to customer needs. For details, please contact the company's customer service or regional sales manager.



Lethen broth swab sampler: pre-filled 4ml or 10ml letheen broth.

Description: In the case of neutralization of disinfectants such as quaternary ammonium is quite important to the antimicrobial activity, it is recommended to use Lethen Broth to test preservative materials.