

## Concentrated Liquid Chromogenic Culture Medium



### Concentrated liquid chromogenic culture medium ST-SSP

The sample solution in ST-25 is discarded by the remaining 5 ml, and then 1.65 ml of the concentrated liquid chromogenic culture medium (ST-SSP) is completely injected to ST-25 for culture.

### Performance Characteristics

ST-SSP can simultaneously detect coliforms and Escherichia coli for qualitative identification.

Name	Description	Packing	Validity period
ST-SSP	Coliforms and Escherichia coli for qualitative identification	10pcs/bag, 300pcs/box	1year

## Operation and Coloring



Add SSP concentrated medium



Primary color – negative



Blue – coliforms positive



Fluorescence -- Escherichia coli positive