

## Sterile Sample Bag



### Sterile Sample Bag

Model	Dimension	Thickness	Capacity	Material	Operating temperature	Packing	Sterilization Method	Validity Period
AK2310	23x10cm	80μm	280ml	Polyethylene	-50~85℃	500 pcs (20 pcs / Ziplock bag x25)	Gamma ray sterilization	3 years
AK2314	23x14cm		550ml					
AK3218	32x18cm		1200ml					

**Precautions:**

1. Do not use for other application than sample sampling, collection, and storage.
2. Keep the sample temperature within the operating temperature range.
3. Used to freeze liquids after which handle and contact carefully.
4. When placing the sample directly in homogenization, be sure to release the sealing strip before homogenizing.
5. After use, remove the contents and dispose of them as specified.
6. In order to improve product quality, specifications and price changes will be notified separately.

**Operation Instruction**

Step 1




Step 2



Step 3



Step 4



Step 5